

MonoPlus®

Absorbable sutures

DESCRIPTION

MonoPlus® is a violet long-term absorbable synthetic monofilament suture made of polydioxanone. MonoPlus® offers good handling combined with great strength. Because of its slow degradation profile, MonoPlus® is ideal for all the indications where an extended wound support of more than 4 weeks is desirable.

INDICATION

- Orthopaedics (tendons, joint capsule)
- Paediatric cardiovascular surgery
- General surgery
- Abdominal wall closure



ADVANTAGES

- High knot tensile strength
- Outstanding strength retention for extended wound support
- Very flexible
- Pliable and easy to knot
- Smooth passage through tissue

MonoPlus®

Product overview

Structure

Monofilament



Colour

Violet

Chemical composition

Polydioxanone

Coating

Uncoated

Origin

Synthetic

Sizes

USP 7/0 (metric 0.5) to USP 2 (metric 5)

Type of absorption

By hydrolysis

Mass absorption

180 up to 210 days

Sterilisation

Ethylene oxide