

# Monomax<sup>®</sup>


## Absorbable Sutures

### DESCRIPTION

Monomax<sup>®</sup> is an extra long-term absorbable monofilament synthetic suture made of poly-4-hydroxybutyrate. Monomax<sup>®</sup> degrades by hydrolysis and enzymatic reaction and provides a predictable and reliable longterm absorption.

### INDICATION

- Closure of the abdominal wall



### ADVANTAGES

- Outstanding degradation profile for abdominal wall closure
- "Shock-Z" – unique extra elasticity
- Unique handling features
- Good adaptation to wound edges

## Monomax<sup>®</sup>

### Product overview

Structure

Monofilament



Colour

Violet

Chemical composition

Poly-4 hydroxybutyrate

Coating

Uncoated

Origin

Synthetic

Sizes

USP 3/0 (metric 2) to USP 1 (metric 4)

Type of absorption

By hydrolysis and enzymatic reaction

Sterilisation

Ethylene oxide