

Monosyn® Quick

Absorbable sutures

DESCRIPTION

Monosyn® Quick is an undyed short-term absorbable synthetic monofilament suture made of glyconate for the approximation of soft tissue with a fast healing process. Monosyn® Quick is recommended for all soft tissue with a fast healing process.

INDICATION

- Plastic surgery
- Gynaecology / Obstetrics
- Urology and Paediatrics
- Odontology and ENT
- Maxillofacial



ADVANTAGES

- Short-term absorption profile

Monosyn® Quick

Product overview

Structure	Monofilament
Colour	Undyed
Chemical composition	Glyconate (72 % glycolide, 14 % trimethylene carbonate, 14 % ϵ -caprolactone)
Coating	Uncoated
Origin	Synthetic
Sizes	USP 6/0 (metric 0.7) to USP 1 (metric 4)
Type of absorption	By hydrolysis
Mass absorption	56 days
Sterilisation	Ethylene oxide