

# Novosyn® QUICK

## Absorbable sutures

### DESCRIPTION

Novosyn® Quick is a polyglactin 910 braided suture with fast degradation profile.

Novosyn® Quick is intended for use in soft tissue approximation when only short-term wound support is required and when the rapid absorption of the suture would be beneficial.

### INDICATION

- Gynaecology / Obstetrics (e. g. episiotomies)
- Skin closure (intracutaneous, subcutaneous)
- Urology
- Odontology (e. g. oral mucosa)
- Paediatric surgery
- Ligatures



### ADVANTAGES

- Fast and predictable degradation profile
- Smooth surface
- High initial tensile strength
- Good knot security
- High flexibility and pliability
- Excellent handling properties
- Quick mass absorption

## Novosyn® QUICK

### Product overview

Structure Braided



Colour Undyed

Chemical composition Polyglactin 910

Coating Polyglactin 370 + calcium stearate

Origin Synthetic

Sizes USP 6/0 (metric 0.7) - USP 2 (metric 5)

Type of absorption in approx. 42 days

Sterilisation Gamma-irradiation