

# Silkam®

## Non-absorbable sutures

### DESCRIPTION

Silkam® is a non-absorbable, braided and coated suture, made of natural silk filaments. Silkam® is basically indicated for applications where it is surgical practice to employ non-absorbable suture material. Silkam® is available in black and white.


### INDICATION

- General surgery
- Skin closure
- Oral surgery
- Neurosurgery
- Ophthalmic surgery
- Ligatures



## Silkam®

### Product overview

Structure	Braided
	
Colour	Black or white
Chemical composition	Natural silk filaments
Coating	Refined paraffin wax
Origin	Natural
Sizes	USP 8/0 (metric 0.4) to USP 6 (metric 8)
Type of absorption	Non-absorbable
Sterilisation	Gamma-irradiation