

Steelex[®]

Non-absorbable sutures

DESCRIPTION

Steelex[®] is a non-absorbable, twisted or monofilament suture, made of corrosion-resistant steel for orthopaedic surgery and sternal closure.

INDICATION

- Orthopaedics
- Sternal Closure



ADVANTAGES

- Exceptional tensile strength
- Excellent tissue compatibility

Steelex[®]

Product overview

Structure Monofilament or twisted



Colour metallic

Chemical composition Stainless steel

Coating Uncoated

Origin Metal

Sizes USP 5/0 (metric 1) to USP 7 (metric 9)

Type of absorption Non-absorbable

Sterilisation Gamma-irradiation