

# Trelon<sup>®</sup>


## Non-absorbable sutures

### DESCRIPTION

Trelon<sup>®</sup> is a non-absorbable synthetic suture made of braided polyamide 6/6 with a silicone coating. It is available in black.

### INDICATION

- General surgery
- Ligatures
- Ophthalmic surgery



## Trelon<sup>®</sup>

### Product overview

Structure Braided



Colour Black

Chemical composition Polyamide

Coating Silicone

Origin Synthetic

Sizes USP 6/0 (metric 0.7) to USP 3 (metric 6)

Type of absorption Non-absorbable

Sterilisation Gamma-irradiation