

DESCRIPTION

In addition to traditional skin suturing, surgeons can use a tissue adhesive.

The tissue adhesive Histoacryl® consists of monomeric n-butyl-2-cyanoacrylate, which polymerises quickly in connection with tissue fluid. Histoacryl® is available in two colors: translucent, special for facial application, and blue, which enables an easy control over the quantity applied.

INDICATIONS

- Closure of Minor Skin Wounds Without Suturing
- Closure of Skin in Endoscopic Incisions

ADVANTAGES

WOUND CLOSURE

- Wound closure in one minute
- Anti-bacterial barrier protects the wound
- Storage at room temperature
- Excellent cosmetic results
- Easy-to-use procedure



Color	Content / box	Cat. No.
Blue	10 ampoules per box	TS1050044-FP
Translucent	10 ampoules per box	TS1050071-FP