

Histoacryl® Flexible

IMPROVED CLOSURE OF SURGICAL INCISIONS

DESCRIPTION

The topical tissue adhesive and microbiological barrier product Histoacryl® Flexible consists of n-butyl-2-cyanoacrylate and a softener. For easy assessment of the thickness of the layer which has been applied, Histoacryl® Flexible is colored blue.

INDICATIONS

- Closure of minimal-tension skin wounds from clean surgical incisions and simple, thoroughly cleansed, trauma-induced lacerations
- Microbiological barrier in the closure of clean surgical skin wounds and simple, thoroughly cleansed, trauma-induced lacerations

ADVANTAGES

- **Flexibility**
Allows closure of incisions up to 25 cm
- **Microbial barrier**
Polymerized Histoacryl® Flexible adhesive films are an effective microbial barrier
- **Ease of use**
Ready to use product storable at temperatures below 25°C
- **Accuracy**
Special tip permits a better control and improved application
- **Excellent cosmetic results**
Cyanoacrylate based topical skin adhesives yield excellent cosmetic results



Color	Content / box	Cat. No.
Blue	10 ampoules per box	1051260-P